



# User Information Sheet

Ronco Item No : **956, 966, 976, 986, 996**

Brand Name : **RONCO Blurite Plus (PF)**

Description : Nitrile Examination Glove  
 Material : Nitrile  
 Donning : Powder Free  
 Characteristics : Disposable, Non-Sterile, Ambidextrous  
Accelerator Free  
 Colour : Dark Blue  
 Texture : Lightly Textured



Specification :

Size :		XS	S	M	L	XL
Item Number :		<b>956</b>	<b>966</b>	<b>976</b>	<b>986</b>	<b>996</b>
Length : (mm)	Minimum	230	230	230	230	230
Width : (mm)	± 10 mm	75	80	95	110	120
Thickness :		4 mil				
Physical Properties :						
	Tensile strength : minimum	11 MPa				
	Elongation : minimum	500%				
Max. powder level:		2 mg/glove				
AQL ( Pinholes )		AQL 1.5				

Packaging : C1M = 100 gloves x 10 dispensers = 1,000 gloves per case

- Product Features :
- Chemo Tested
  - Accelerator-Free
  - Latex Free
  - Superior elasticity
  - Ease of finger movement
  - Reduces hand fatigue
  - Resists a wide variety of chemicals
  - Excellent strength
  - Textured
  - Beaded cuffs

- Applications :
- Food processing
  - Food service
  - Hospitality
  - Meat packing
  - Dental
  - Medical
  - Pharmaceutical
  - Veterinary
  - Automotive
  - Beauty / Aesthetics
  - Electronics
  - Environmental
  - Industrial
  - Laboratory
  - Printing

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F ( 37°C )

Special Notes : Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).